

**1. Technische Daten**

Ausgabekapazität bei einer Zapfrate von 4 Getränken a 0,3 Liter pro Minute**	316
Eisbankgewicht in kg	18
Eisbankkapazität in kcal	1440
Eisaufbau in Minuten ohne Python	200
Anschlußspannung	220V/60Hz
Leistungsaufnahme in Watt	1050
Kompressorleistung in Watt*	509
Kältemittel R 134 a in kg	0,350
Karbonatorpumpenleistung in Liter/Stunde bei 10 bar	280
Umwälzpumpenleistung in Liter/Stunde bei 2 bar	320
Kühlleistung/Eisbankleistung in Watt	517
in kcal/h	445
Anzahl Kühlschlangen	
Sirup	6
Premix	1
Stillwasser	1
Abmessungen in mm	
Höhe	595
Breite	780
Tiefe	433
Versandgewicht in kg	85

**1. Technical Data**

Output capacity at a tap rate of 4 drinks of 0.3 liters each per minute**
Weight of ice bank in kg
Ice bank performance in kcal
Ice build up in minutes without python
Supply voltage
Power consumption in watt
Compressor output in watt*
Refrigerant R134 a in kg
Carbonator pump output in liter hour at 10 bar
Circulation pump output in liter/hour at 2 bar
Cooling / ice bank performance in watt
in kcal/h
Number of cooling coils
Syrup
Premix
Still water
Dimensions in mm
Height
Width
Depth
Shipping weight in kg

\* bei -10°C Verdampfungstemperatur  
\*\*mit 10 m SC Python

\* at -10°C evaporation temperature  
\*\* with 10 m SC python

Kühlleistungen und Ausgabekapazität bei 24°C Umgebungstemperaturen und Wasser- bzw. Sirupeingangstemperaturen von 24°C und Getränkeausgangstemperaturen von unter 5°C. Bei Verwendung von Cornelius Pythons muß pro lfd. Meter mit einem Kühlleistungsverlust von 13 kcal/h gerechnet werden. Änderungen vorbehalten.

Cooling capacities and output capacity at 24°C ambient temperatures and water or syrup inlet temperatures of 24°C and beverage outlet temperatures of less than 5°C. When Cornelius pythons are used, a cooling loss of 13 kcal/hour per running meter must be included in calculations. We reserve the right to make modifications.

**2. Position der Eisbankelektrode**

**2. Position of the Icebank Probe**